

Let's Play Sport Wheelchair Measurement Guide

3 Easy Steps to Determine appropriate size:

- Measure Seat Width and Depth - Determine seat width by measuring the person's hips from one side to the other, in a straight line. Determine seat depth by measuring from the back of the hip to the back of the knee of the person while seated. Then, subtract one inch from this measurement to choose the correct seat depth.
- Footrest - To determine the length of the footrest, measure the distance from the back of the person's knee to the heel of their foot.
- Back Height - Measure the distance from the person's armpit down to the seat.

